



FOR USE IN THE CATH LAB, EP, MRI AND
RADIOLOGY OR WHEREVER IMAGING OCCURS



FEATURES

- Radiolucent Ag/AgCL sensor
- Hypoallergenic adhesive
- Latex Free
- Available in economically packaged, no-waste quantities

1.800.669.6905 | Vermed.com

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

MR Conditional



Latex-free



Wet Gel



Solid Gel





ClearScan™



Radiolucent Electrodes for ECG Monitoring

Keep your imaging free of obstruction by using Vermed ClearScan™ Electrodes. Designed with a carbon-filled plastic snap that will not interfere with your Nuclear Imaging Systems.



A10004-5 - 1.5" / 38.1MM - ROUND FOAM ELECTRODE WITH CARBON SNAP

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	MR Conditional
A10004-5	5	5	100	600		
A10004-60	5	60	n/a	600		



A10004-5T - 1.625" / 41.275MM - ROUND TAPE ELECTRODE WITH CARBON SNAP

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	MR Conditional
A10004-5T	5	5	100	600		
A10004-60T	5	60	n/a	600		

A10006-5RT - 1.875" / 47.625MM - ROUND FOAM ELECTRODE WITH CARBON SNAP

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	MR Conditional
A10006-5RT	5	5	100	600		
A1006-1-60RT	1	60	n/a	600		


A10025-60 - 1.875" / 47.625MM - ROUND FOAM ELECTRODE WITH CARBON SNAP

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	MR Conditional
A10025-60	5	60	n/a	600		

A10031-60RT - 1.875" / 47.625MM - ROUND FOAM ELECTRODE WITH CARBON SNAP

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10031-60RT	5	60	n/a	600		

A10032-60RT - 1.875" / 47.625MM - ROUND TAPE ELECTRODE WITH CARBON SNAP

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10032-60RT	5	60	n/a	600		

A10004-5

Foam | Wet Gel

A10031-60RT

Foam | Solid Gel

A10004-5T

Tape | Wet Gel

A10032-60RT

Tape | Solid Gel

A10006-5RT / A10025-60

Foam | Wet Gel

