

FEATURES

- Labeled specifically for EMS
- Packaged in convenient “no waste” quantities for EMS
- 2.25” (57.15 mm) diameter provides large adhesive area
- High chloride wet gel provides fast, clear traces
- Can be defibrillated
- Hypoallergenic and latex free

1.800.669.6905 | Vermed.com

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

Latex-free



Wet Gel




EMS Electrodes

When every second counts

Vermed EMS Electrodes are designed and packaged specifically for Emergency Services. Large size and aggressive adhesive ensure quick placement and fast, accurate readings. All Vermed Electrodes use the highest quality materials and meet or exceed AAMI specifications for disposable ECG Electrodes


A10003 - 2.25" / 57.15MM - ROUND FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10003-3	3	3	75	600	
A10003-4-60	4	60	n/a	600	
A10003-4-48	4	48	n/a	48	


A10084-4 - 1.625" / 41.275MM - ROUND TAPE ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10084-4	4	4	100	600	

A10005-4NB - 1.5" / 38.1MM - ROUND FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10005-4NB	4	4	n/a	600	

A10012-6S - 1.25" X 1.875" / 31.75MM X 47.625MM - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10012-6S	6	6	n/a	600	

A10003
Foam | Wet Gel



A10084-4
Tape | Wet Gel



A10005-4NB
Foam | Wet Gel



A10012-6S
Foam | Wet Gel

