

12 LEAD ECG KIT

V1 - V2:

Small rectangular electrode with two sensors

V3 - V6:


Small teardrop shaped electrode for placement in tight areas

RA, LA, RR, LL:

Large round electrodes for superior adhesion

SINGLE USE PACKS

Packaged in no-waste quantities



CASE

FDA / GS1 Compliant Labels



1.800.669.6905 | Vermed.com

Patient Monitoring & Medical Surgical Consumables

Latex-free



Wet Gel




PerformancePlus™ Electrodes

Stress Test Kit


Our unique single use pouch containing 10 electrodes is perfect for diagnostic stress testing. Electrodes come in varying sizes for ease of placement in designated positions.

A10015-10S - STRESS TEST KIT - FOR 10 LEAD PLACEMENT

10 ELECTRODES/ KIT | 50 KITS / CASE

Size	Shape	Per Pouch	Wet Gel
2.25" / 57.15mm	Round	4	
1.5" x 38.1mm	Teardrop	4	
1.75" x 1.625" / 44.43mm x 41.275mm	Square, dual element	1	

WHAT'S IN A KIT?



4 2.25" / 57.15mm
Foam | Wet Gel

4 1.5" / 38.1mm Teardrop
Foam | Wet Gel

1 1.75" x 1.625" /
44.45mm x 41.275mm
Perforated | Foam | Wet Gel



Perforation allows separation when desired