




FEATURES

- Ag/AgCL sensor
- High chloride wet gel hydrates the skin and produces fast and clear traces
- Aggressive adhesive keeps electrodes on patient - no motion artifact or repeat testing
- Hypoallergenic adhesive



SINGLE USE PACKS

Packaged in no-waste quantities



INNER SLEEVES

Single use packs are packaged in inner sleeves within the case

RESEALABLE BULK POUCH

Contains 30 or more electrodes



CASE
FDA / GS1 Compliant Labels



1.800.669.6905 | Vermed.com

Patient Monitoring & Medical Surgical Consumables





PerformancePlus™ Electrodes

Electrodes for Stress and Holter Testing

In critical diagnostic situations, choose Vermed PerformancePlus™ Electrodes to ensure quality readings, limited artifacting, and no repeat procedures.

A10005 - 1.5"/38.1mm - TEARDROP FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10005	5	5	100	600		 24 Hours
A10005NB	5	5	n/a	600		
A10005-1-60	1	60	n/a	600		
A10005-4NB	4	4	n/a	600		
A10005-7NB	7	7	n/a	700		
A10005-60	5	60	n/a	600		



A10006-5 - 1.875"/47.625mm - TEARDROP FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10006-5	5	5	100	600		 72 Hours
A10006-5-50	5	50	n/a	600		
A10006-1-60S	1	60	n/a	600		



A10006-5T - 2.25"/57.15mm - ROUND TAPE ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10006-5T	5	5	100	600		 48 Hours
A1006-1-60TS	1	60	n/a	600		





A10007-5 - 2"/50.8mm - ROUND CLOTH ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10007-5	5	5	100	600		 72 Hours
A10007-3	3	3	n/a	600		





A10012-5S - 1.875" x 1.25" / 47.625mm x 31.75mm - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10012-3S-60	3	3	n/a	30		 48 Hours
A10012-5S	5	5	100	600		
A10012-4-60S	4	60	n/a	600		
A10012-5S	5	5	100	600		
A10012-6S	6	6	n/a	600		
A10012-6S-60	6	6	n/a	60		
A10012-60S	5	60	n/a	600		
A10012-10S	5	10	60	600		



A10101-3 - 1.875" x 1.25" / 47.625mm x 31.75mm - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10101-3	3	3	108	1080		 24 Hours
A10101-10	5	10	100	1000		

