



Perforation allows separation
when desired



FEATURES

- Dual Element
- Silver / Silver Chloride universal snaps
- Convenient small size - 1.625" x 3.25" (41.275 x 82.55mm)
- Non-invasive
- Pre-gelled and disposable
- 1/16" Aggressive foam

1.800.669.6905 | Vermed.com

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

Latex-free



Wet Gel



PNS Dual Element™

Reliable Testing of the effectiveness of nerve blockers

For accurate and reliable testing of the effectiveness of nerve blockers, Vermed PNS Dual Element™ electrodes are a must. They are pre-gelled, disposable, and feature universal Silver / Silver Chloride snaps to provide positive stimulation. Most importantly, they are made with the same care and high standards for which all Vermed products are known.

A10011 - 3.25" X 1.625" / 82.55MM X 41.275MM - RECTANGULAR DUAL ELEMENT ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case
A10011	1	1	50	300

A10011
Foam | Wet Gel

