



## FEATURES

- High Chloride wet gel provides fast, clear traces
- Self-adhering electrodes help reduce muscle artifact or repeat testing
- Available in economical, pliable, easy-to-peel vinyl tape or easy placement, more adhesive, occlusive foam
- Can be defibrillated
- Hypoallergenic and Latex Free

**1.800.669.6905 | Vermed.com**

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

Latex-free



Wet Gel



Solid Gel



# SilveRest™ Electrodes

## Adult Diagnostic Electrodes for Resting EKG and Echo Sound

Vermed SilveRest™ Electrodes are the standard in diagnostic electrodes. Superior Silver / Silver Chloride technology and quality construction produce an electrode that performs every time. Excellent adhesion and clear traces mean no repeat testing.


### A10009 - 1.375" / 34.925MM - ROUND VINYL TAPE ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10009	10	10	200	1200	
A10009-100	10	100	n/a	1000	

### A10009-100F - 1.375" / 34.925MM - ROUND FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10009-100F	10	100	n/a	1000	

### A10013-100S - 1" X 1.5" / 25.4MM X 38.1MM - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10013-100S	5	100	n/a	1000	

### A10049 - 0.8125" X 1.75" / 20.638MM X 44.45MM - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10049	5	50	500	2000	

• Quality Construction | Superior Performance | Economical Testing •

#### A10009

Vinyl | Wet Gel




#### A10009F

Foam | Wet Gel



#### A10013

Foam | Wet Gel



#### A10049

Foam | Solid Gel

