

FEATURES

Pre-Wired Electrodes:

We offer a complete line of pre-wired electrodes for use in PICU/NICU and other neonatal applications. Available in various sizes as well as metallic or radiolucent lead wires.

Probe Covers:

Our thermal probe covers use a delicate adhesive designed for tender infant skin. The highly reflective surface deflects heat from radiant warmers.

1.800.669.6905 | Vermed.com

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

MR Conditional Latex-free Solid Gel



TenderTrode™



Pre-Wired Electrodes & Thermal Probe Covers

Vermed offers a complete line of pre-wired electrodes for use in the PICU/ NICU and other neonatal applications. Available in various sizes as well as metallic or radiolucent lead wires.


A10057-S - 0.875" / 22.225MM - SQUARE PRE-WIRED FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10057-S	3	3	n/a	300	



A10061-SRT - 0.875" / 22.225MM - SQUARE PREWIRED CLOTH ELECTRODE - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10061-SRT	3	3	n/a	300		


A10010-S - 0.875" X 1.875" / 22.225MM X 47.625MM - RECT. PREWIRED CLOTH ELECTRODE


Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10010-S	3	3	n/a	300	

A10060-SRT - 0.875" X 1.875" / 22.225MM X 47.625MM - RECT. PREWIRED CLOTH ELECTRODE - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	MR Conditional
A10060-SRT	3	3	n/a	300		

A10008 - THERMAL PROBE COVER

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10008	1	50	n/a	500	



A10057-S
Foam | Solid Gel
A10061-SRT
Radiolucent
Cloth | Solid Gel

A10010-S
Cloth | Solid Gel
A10060-SRT
Radiolucent
Cloth | Solid Gel

A10008
Cloth | Foil

