




FEATURES

5 DAY WEAR TIME, GUARANTEED!

- Hydrocolloid is gentle on delicate skin
- Provides long-term attachment with no skin irritation
- Resists breakdown from moisture or heat
- Leaves no residue on skin
- Save nursing time and resources with fewer electrode changes
- Safe on neonatal skin
- Designed specifically for NICU
- Radiolucent lead wires



A10055-SRT (3 Lead)
A10077-SRT (5 Lead)
Cloth | Solid Gel


A10076
Hydrocolloid Temperature
Probe Cover

A10058-SRT
Cloth | Solid Gel

A10086
Mini Cannula Holder

PRE-WIRED ELECTRODE ANATOMY


Each electrode is designed with gentle, hypoallergenic adhesives and soft, conforming materials to ensure comfort and effectiveness on young patients.



**Enlarged for display*

SINGLE USE PACKS

Packaged in no-waste quantities



CASE
FDA / GS1 Compliant Labels



1.800.669.6905 | Vermed.com

Patient Monitoring & Medical Surgical Consumables

Latex-free



Solid Gel





TenderTrode Plus™



Pre-Wired Electrodes & Thermal Probe Covers

Vermed's TenderTrode Plus™ line of hydrocolloid products are designed for the unique needs of the neonatal market. Each product uses a hydrocolloid instead of aggressive tape adhesives to prevent irritation to tender skin while providing long-term attachment. This, coupled with the use of soft conforming materials, ensures patient comfort and product effectiveness in the NICU.



A10055-SRT - 1" / 25.4MM - ROUND PRE-WIRED CLOTH ELECTRODE - 3 LEAD - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10055-SRT	3	3	n/a	60		 120 Hours



A10077-SRT - 1" / 25.4MM - ROUND PRE-WIRED CLOTH ELECTRODE - 5 LEAD - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10077-SRT	5	5	n/a	100		 120 Hours



A10058-SRT - 1" X 0.375" / 25.4MM X 9.525MM - RECT. PRE-WIRED CLOTH ELECTRODE - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10058-SRT	3	3	n/a	60		 120 Hours



A10076 - 1" / 25.4MM - HYDROCOLLOID TEMPERATURE PROBE COVER

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10076	1	50	n/a	500		 120 Hours

A10086 - MINI CANNULA HOLDER FOR NEONATES

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10086	1	n/a	n/a	100		 120 Hours

A10085 - CANNULA HOLDER

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10085	1	n/a	n/a	100		 120 Hours

Clinical Feedback

Recent clinical feedback indicates a **98% reduction in faulty incubator and monitor alarms** tripped when standard electrodes peel off in high heat and humidity.

TenderTrode™ Plus is now **leading the way** in neonatal patient care. **Significant cost savings** are evident due to limited replacement, a function of our patented hydrocolloid adhesive border. These electrodes have a 5 day wear time, guaranteed!

Reducing the removal and replacement of standard monitoring electrodes **ensures the body is focused on essential growth and repair.** Reduced monitor alarms ensures the comfort of patients and parents alike while reducing nursing staff workload.