




**Compares to both the Ambu<sup>®</sup> BlueSensor and Skintact<sup>®</sup> Tab 'n' Snap.**

## SimplySnap

Offset Electrodes

- Offset connector minimizes artifact from patient & cable movements
- Clinicians can attach cables while electrodes are on patient without compromising gel well
- Wet gel provides quick reading & stabilization
- Intuitive application process



*Stainless Steel, Carbon, Offset Tab, and pediatric styles available.*

## VersaTrode<sup>+</sup>

Snap 'n Tab Electrodes

- Electrode can be used with either a snap or an alligator clip
- MR Conditional / Radiolucent snap
- Suitable for all diagnostic procedures
- Stays on patient throughout multiple departments



**1.800.669.6905 | Vermed.com**

Patient Monitoring & Medical Surgical Consumables

Latex-free





MR Conditional



Wet Gel





## SimplySnap

Offset Electrodes

### PACKAGING



#### RESEALABLE BULK POUCH

Resealable pouch keeps electrodes fresh after opening



#### CASE

FDA / GS1 Compliant Labels

							<b>CROSS REF</b>	
Description	Product #	Reference	Size	Snap	Per Pouch	Per Box	SkinTact	Ambu
SimplySnap Offset Electrode	4010425C	OTL-5235-S	52 x 35	Stainless	60	1200	FS-VB01	Ambu R-00-S Ambu T-00-S Ambu SP-00-S
	32029077	OTL-5235-S-5	52 x 35	Stainless	5	1000		
SimplySnap Radiolucent Offset Electrode	4010452C	OTL-5235-SRT	52 x 35	Carbon	60	1200	N/A	N/A
	32029078	OTL-5235-SRT-5	52 x 35	Carbon	5	1000		
Offset Electrode No Snap	4010423C	OTL-5235	52 x 35	None	100	1500	FS-VB00	Ambu SU-00-C
	32029079	OTL-5235-5	52 x 35	None	5	1000		
SimplySnap Pediatric Offset Electrode	4010453C	OTL-4227-S	42 x 27	Stainless	60	1200	F-VK01	Ambu P-00-S/50
	32029080	OTL-4227-S-5	42 x 27	Stainless	5	1000		
VersaTrode+ Snap 'n Tab	4010429C	TSL-4532	45 x 32	Carbon	60	1200	FS-WA10C/5	Ambu BS-Tab-2300/100
	32029076	TSL-4532-5	45 x 32	Carbon	5	1000		

The trademarks identified herein are the trademarks of the respective OEMs, and not Vermed. Vermed disclaims any affiliation, connection or association between its products and those of the respective OEMs.